

PrimePad® Series Product Matrix



Code Prefixes:
PrimePad-(code suffix)
PrimePadPlus-(code suffix)
 Example: PrimePadPlus-K2

Code Suffix	Desc.	Case Pack	LxW	Thickness	Table Rail Straps	Body Strap Set	Draw Sheet	Head Cradle	Arm Protection Options		
									Sleeves	Arm Protectors	Wings
1	Trendelenburg Pad	10	112 cm x 58 cm	2.54 cm							
1S	Trendelenburg Pad	10	112 cm x 58 cm	2.54 cm	●						
2	Trendelenburg Pad	4	112 cm x 58 cm	5.08 cm							
2S	Trendelenburg Pad	4	112 cm x 58 cm	5.08 cm	●						
3	Trendelenburg Pad	4	112 cm x 58 cm	2.54 cm to 5.08 cm							
3S	Trendelenburg Pad	4	112 cm x 58 cm	2.54 cm to 5.08 cm	●						
BBO	Bean Bag Overlay	6	101.6 cm x 96.52 cm	2.54 cm							
CVP	Cardiovascular Pad	4	183 cm x 58 cm	2.54 cm							
K1	Trendelenburg Pad	10	112 cm x 58 cm	2.54 cm	●	●					
K2	Trendelenburg Pad	10	112 cm x 58 cm	2.54 cm	●	●	●				
K3	Trendelenburg Pad	3	112 cm x 58 cm	2.54 cm	●	●			●		
K6	Trendelenburg Pad	4	112 cm x 58 cm	2.54 cm	●	●	●		●		
K7	Trendelenburg Pad	3	112 cm x 58 cm	2.54 cm	●	●	●	●	●		
K8	Trendelenburg Pad w/integrated arm tucking/protection	3	112 cm x 58 cm	2.54 cm	●	●	●	●			61 cm x 41 cm
K9	Trendelenburg Pad	10	112 cm x 58 cm	2.54 cm	●		●				
K10	Trendelenburg Pad	4	112 cm x 58 cm	5.08 cm	●	●	●				
K12	Trendelenburg Pad	3	112 cm x 58 cm	2.54 cm	●	●	●	●			
K14	Trendelenburg Pad	8	112 cm x 58 cm	2.54 cm	●	●	●			●	
K15	Trendelenburg Pad	3	112 cm x 58 cm	2.54 cm	●	●	●	●		●	
K16	Trendelenburg Pad w/integrated arm tucking/protection	4	112 cm x 58 cm	2.54 cm	●	●					61 cm x 41 cm
W1	Trendelenburg Pad w/integrated arm tucking/protection	4	112 cm x 58 cm	2.54 cm	●						61 cm x 61 cm
W2	Trendelenburg Pad w/integrated arm tucking/protection	4	112 cm x 58 cm	2.54 cm	●						61 cm x 41 cm
PED	Pediatric Trendelenburg Pad	10	61 cm x 51 cm	2.54 cm							
SQ1	Square Positioning Foam	25	30.5 cm x 30.5 cm	2.54 cm							