

PrimePad® Series Product Matrix



Code Prefixes:
PrimePad-(code suffix)
PrimePadPlus-(code suffix)
 Example: PrimePadPlus-K2

Code Suffix	Desc.	Case Pack	LxW	Thickness	Table Rail Straps	Body Strap Set	Draw Sheet	Head Cradle	Arm Protection Options		
									Sleeves	Arm Protectors	Wings
1	Trendelenburg Pad	10	44"x23"	1"							
1S	Trendelenburg Pad	10	44"x23"	1"	●						
2	Trendelenburg Pad	4	44"x23"	2"							
2S	Trendelenburg Pad	4	44"x23"	2"	●						
3	Trendelenburg Pad	4	44"x23"	1" to 2"							
3S	Trendelenburg Pad	4	44"x23"	1" to 2"	●						
BBO	Bean Bag Overlay	6	40"x38"	1"							
CVP	Cardiovascular Pad	4	72"x23"	1"							
K1	Trendelenburg Pad	10	44"x23"	1"	●	●					
K2	Trendelenburg Pad	10	44"x23"	1"	●	●	●				
K3	Trendelenburg Pad	3	44"x23"	1"	●	●			●		
K6	Trendelenburg Pad	4	44"x23"	1"	●	●	●		●		
K7	Trendelenburg Pad	3	44"x23"	1"	●	●	●	●	●		
K8	Trendelenburg Pad w/integrated arm tucking/protection	3	44"x23"	1"	●	●	●	●			24" x 16"
K9	Trendelenburg Pad	10	44"x23"	1"	●		●				
K10	Trendelenburg Pad	4	44"x23"	2"	●	●	●				
K12	Trendelenburg Pad	3	44"x23"	1"	●	●	●	●			
K14	Trendelenburg Pad	8	44"x23"	1"	●	●	●			●	
K15	Trendelenburg Pad	3	44"x23"	1"	●	●	●	●		●	
K16	Trendelenburg Pad w/integrated arm tucking/protection	4	44"x23"	1"	●	●					24"x16"
W1	Trendelenburg Pad w/integrated arm tucking/protection	4	44"x23"	1"	●						24"x24"
W2	Trendelenburg Pad w/integrated arm tucking/protection	4	44"x23"	1"	●						24"x16"
PED	Pediatric Trendelenburg Pad	10	24"x20"	1"							
SQ1	Square Positioning Foam	25	12"x12"	1"							